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Crimping tool T.802


Crimping tool T.802
6971


Warning: We recommend to operate this complete check up of the crimping tool at least every 500 crimping operations.

1. **TIGHTENING OF SCREWS**

Check that all screws are well in place and blocked (see pictures **A** et **B**).

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
DO NOT TRY TO ADJUST THE SCREW shown in picture D, simply ensure it is tightened.

Picture D

2. POSITIONING OF THE CRIMPING SPRING

Check the correct position of the crimping spring # 6598 everywhere well inside the groove: see picture **E**.

Picture E

Spring # 6598

Check that the spring #6598 is not damaged. If damaged, replace it by new spring.

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3. CLEANING OF THE CRIMPING TOOL

In order to guarantee optimum operation: clean the crimping tool at regular intervals by removing any dirt or particulate for instance with a soft tissue moistened with alcohol.

4. GREASING OF THE CRIMPING TOOL

Use Food Grade silicon lubricant to grease the following parts:

- The spring: remove it, coat it with a fine layer of sprayed lubricant and put it back into the crimping tool.
- The axles.
- The inner side of the crimping head.

Food Grade lubricant recommended: KLUBER UNISILKON M2000 in spray or equivalent.

5. STERILIZATION

The sterilization process must not expose the tool to temperatures above 125°C and the use of oxydant products may damage the tool or create pitting.

Cleaning, drying and lubricating the crimping tool are necessary operations in order to keep the tool functioning properly.

IMPORTANT: In case of doubt about performances of the sealing tool, it is recommended to contact our commercial department for assistance.

NOTE: *The technical characteristics of this crimping tool and the information of the user's guide are provided as general information and can be updated without prior notice. TOURNAIRE S.A shall not be liable for any incorrect use of its containers and accessories or for any wrong understanding of the present "maintenance form". In any case, the filler remains the sole responsible to validate its filling and closing process, to operate according to Good Manufacturing Practices and to carry out appropriate controls of the finished product.*

The present maintenance form is available in e-format; contact TOURNAIRE S.A for more information

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